
LXDeviceAPI

LXDeviceAPI Developer Manual

Doc. ID. LXE64 V1

Release Date. 2017-04-25 .

Abstract –API for bio signal measuring devices



Site for LXDeviceAPI : <http://laxtha.net/lxdeviceapi/>



목차

OVERVIEW LXDEVICEAPI	5
FEATURES LXDEVICEAPI	5
GETTING STARTED	6
STEP 1. DOWNLOAD / UNBLOCK / UNZIP LXDEVICEAPI	6
step1.1 Download	6
step1.2 Unblock the zip file	6
step 1.3 unzip	7
STEP 2 . COPY FILES TO YOUR PROJECT FOLDER	7
STEP 3.IMPORTING DLL	7
API STREAM DATA	8
STRUCT : ST_STREAMDATA_LXDAPI	8
USAGES	8
MEMORY MAP	9
Wave_StreamData_CS Data Allocation	9
Event_StreamData_CS Data Allocation	9
ARRAY INDEXING FOR WAVE DATA	10
ARRAY INDEXING FOR EVENT DATA	10
INDEXING EXAMPLE	11
API MESSAGE	12
MESSAGES BY LXDEVICEAPI	12
MESSAGE PARAMETER - WPARAM	12
MESSAGE TYPE VS. MEANING	12
CODE EXAMPLES	13
Code Example 1 . SetMessageDevice	13
Code Example 2. StartStream	13
Code Example 3. Message Handler	13
API FUNCTIONS	14
OPENAPI	15
PARAMETER	15
Return Values	15
Code Example	15
CLOSEAPI	16
Return Values	16
Code Example	16
OPENDevice	17
Parameter	17
Return Values	17
Code Example1	18
Code Example2	18
CLOSEDevice	19
Parameter	19




Return Values	19
Code Example	19
STARTSTREAM	20
Parameter	20
Return Values	20
Code Example	20
STOPTSTREAM	21
Parameter	21
Return Values	21
Code Example	21
GETSTREAMDATA	22
Parameter	22
Return Values	22
Code Example	22
EVENTMARKINGONSTREAM	24
Parameter	24
Return Values	24
Code Example 1	24
Code Example 2	25
SETMESSAGEDEVICE	26
Parameter	26
Return Values	26
Code Example	26
SETSAMPLEFREQUENCY	27
Parameter	27
Return Values	27
Code Example	27
SETDEVICECONTROLPANEL	28
Parameter	28
Return Values	28
Code Example	28
GETFILTERFREQUENCY	29
Parameter	29
Return Values	29
Code Example	30
GETSAMPLEFREQUENCY	31
Parameter	31
Return Values	31
Code Example	31
GETEEGREFELECTRODE	32
Parameter	32
Return Values	32
Code Example	32
CHECKFORUPDATE	33
Parameter	33
Return Values	33
Code Example	33
APPENDIX 1. SUPPORTING DEVICES	34
REVISION HISTORY	35





Overview LXDeviceAPI

The devices can be communicated with your application program via LXDeviceAPI. You can make your own application programs for communicating with devices. LXDeviceAPI provides the standard C libraries and messages.



System Architecture

Features LXDeviceAPI.

- Real time bio signal data streaming from measuring device to your application program.
- Extremely stable and reliable thread for real time streaming data.
- Setting the device parameters (Sampling Frequency, Device Mode Change)
- Providing the self contained UI.
 - Setting the signal filter(LPF, HPF, notch) via “device control panel”.
 - Electrode-Skin Impedance Monitoring.
 - Auto Calibrating the device.
 - Auto Self Updating API.
 - Auto Device Firmware Update.
 - Saving the configuration information into Device.
- API type : DLL (Dynamic Link Library).
- API Making tool : Visual C++ 2015. MFC Regular DLL project.
- Supporting platform : both 32bit and 64bit.
- Supporting OS : Windows 10, 8.1, 8, 7.

Getting Started

Step 1. Download / Unblock / Unzip LXDeviceAPI

step1.1 Download.

Click below link,

LXDeviceAPI_32bit.zip for 32 bit Applications.


https://github.com/LAXTHA/LXDeviceAPI/raw/master/LXDeviceAPI_32bit.zip

LXDeviceAPI_64bit.zip for 64bit Applications.

https://github.com/LAXTHA/LXDeviceAPI/raw/master/LXDeviceAPI_64bit.zip


step1.2 Unblock the zip file..

Before unzip, unblock the zip file should be taken. Click the downloaded zip file by right mouse and click the Properties (red box following left image). Check Unblock and then click OK.



step 1.3 unzip.

After unzip, you can see the files and folder as following.



Step 2 . Copy files to your project folder

Files	Copy to
<ul style="list-style-type: none"> • LXDeviceAPI_Utility (폴더) • LXDeviceAPI.DLL • chartdir60.dll 	copy to the same folder which contains your application exe file.
<ul style="list-style-type: none"> • LXDeviceAPI.LIB • LXDeviceAPI.h 	copy to your application project's source folder.

Step 3.Importing DLL

In your Visual C++ source file, add the codes for implicit linking the library and include the LXDeviceAPI.h

```
#pragma comment(lib,"LXDeviceAPI.lib") // implicit linking.
#include "LXDeviceAPI.h" // for using LXDeviceAPI
```

That's all, now it's ready to use all functions and message from LXDeviceAPI.



API Stream Data

Struct : ST_STREAMDATA_LXDAPI

LXDeviceAPI.h provides the typedef of ST_STREAMDATA_LXDAPI.

```
// struct for stream.
typedef struct _stStreamData_LXDAPI
{
...
    unsigned int*      Event_StreamData_CS; // array for event stream
...
    double*            Wave_StreamData_CS; //array for measured multi channel bio signal
...
}ST_STREAMDATA_LXDAPI;
```

The most important member variables are Wave_StreamData_CS for multi channel bio signal and Event_StreamData_CS for event marking.

double * Wave_StreamData_CS

unsigned int* Event_StreamData_CS

Usages.

Declare the variable for stream data type and feed into OpenDevice as a address of variable.

```
ST_STREAMDATA_LXDAPI stStreamData; // .
int NumSampleReturn = 32 ; // Available Value ㄱ : 1 ~ 128
OpenDevice_LXDeviceAPI( , &stStreamData, NumSampleReturn, ); //
```


LXDeviceAPI dynamic allocates the memory for arrays when you call the OpenDevice. The size of memory is determined by the OpenDevice's parameter int numsample_return which is the number of samples per each message received. LXDeviceAPI deletes the memory if CloseDevice is called.



Memory Map.

Wave_StreamData_CS Data Allocation.


Wave Stream Data. Memory Map & Array Indexing.



- N : Number of Samples
- M : Number of Channels
- **sample_index** : index for samples. from 0 to N-1
- **channel_index** : index for channels. from 0 to M-1

Event_StreamData_CS Data Allocation

Event Stream Data. Memory Map & Array Indexing.



- N : Number of Samples
- **sample_index** : index for samples. from 0 to N-1



Array Indexing for Wave data

The following code shows how to get the one sample for specific channel index and sample index,

```
double one_sample_wave = stStreamData.Wave_StreamData_CS[sample_index + NumSampleReturn *  
channel_index];
```

where,

channel_index : Available value range from 0 to (number of channel – 1). The number of channel is specific to device. In the case of QEEG-64FX, the number of channel is 67 (EEG 64ch + Bipolar 3ch).

sample_index : Available value range from 0 to (NumSampleReturn – 1). The NumSampleReturn is determined by the parameter int numsample_return when you call the OpenDevice.

Array Indexing for Event Data

```
unsigned int one_sample_event = stStreamData.Event_StreamData_CS[sample_index];
```


where,

sample_index : Available value range from 0 to (NumSampleReturn- 1).The NumSampleReturn is determined by the parameter int numsample_return when you call the OpenDevice.



Indexing Example

Array indexing example for two channel signals and event. The chart plots the both signals and event simultaneously.





API Message

Messages by LXDeviceAPI

LXDeviceAPI send message to application program using win32 function SendMessage(...,wParam,lParam) . The message parameter wParam has informations to recognize the message type and what processing should be taken by application program.

Message Parameter - wParam

LXDeviceAPI send a wParam as an unsigned int type(byte size 4). Application program can take the information at each byte as follows; wParam is divided into 4 bytes. Byte0 is lowest byte and Byte3 means highest byte.

Parameter	Written Data	Description.
wParam_Byte3	Device Handling ID	<ul style="list-style-type: none">Device Handling ID = wParam_Byte3*256 + wParam_Byte2. range : 1~65535The message is sent by device correspond to "Device Handling ID".Related : Device Handling ID can be retrived by return value of OpenDevice.
wParam_Byte2		
wParam_Byte1	Message Type ID.	<ul style="list-style-type: none">Identificating the message.available values : MSGTYPEID#_DEVICE_LXDAPI, MSGTYPEID#_API_LXDAPI predefined on LXDeviceAPI.h.Related Function : SetMessageDevice
wParam_Byte0	Message Type Sub ID.	<ul style="list-style-type: none">sub message type under Message Type ID.

Message Type vs. Meaning

wParam Byte 1 (Message Type ID)	wParam Byte 0 (Message Type Sub ID)	Description
MSGTYPEID0_DEVICE_LXDAPI	0	The number of sampling data is ready. App's message handler must call the function GetStreamData to get the stream data from LXDeviceAPI. The "number of sampling" is defined by int numsample return when calling the function OpenDevice.



Code Examples

Code Example 1 . SetMessageDevice

```
#define WM_STREAM_DEVICE WM_USER+203 // Define Message to get message from LXDeviceAPI.

void CLXDeviceAPI_Sample1View::OnMenuSetmessagedevice()
{
    SetMessageDevice_LXDeviceAPI(m_iDeviceHandlingID, MSGTYPEID0_DEVICE_LXDAPI ,this->m_hwnd,
WM_STREAM_DEVICE,1); //
}
```

Code Example 2. StartStream

```
void CLXDeviceAPI_Sample1View::OnMenuStartstream()
{
    StartStream_LXDeviceAPI(m_iDeviceHandlingID);
}
```

Code Example 3. Message Handler

After calling the function StartStream, the LXDeviceAPI sends a message when the number of sample data is ready. Your app's message handler must call the function GetStreamData to get the stream data from LXDeviceAPI. The "number of sampling" is defined by int numsample_return when calling the function OpenDevice. If you stop the message from LXDeviceAPI, call the function StopStream or call the SetMessageDevice as last parameter = 0.

```
/* Real Time Acquisition */
afx_msg LRESULT CLXDeviceAPI_Sample1View::OnStreamData(WPARAM wParam, LPARAM lParam)
{
    unsigned int uintWPARAM = (unsigned int)wParam;


    unsigned char msgtype_id = (unsigned char)(uintWPARAM >> 8); //get the lowest 2'nd byte(message type id).
    unsigned char msgtype_subid = (unsigned char)(uintWPARAM); //get the lowest 1st byte(message type sub id).

    switch (msgtype_id)
    {
        case MSGTYPEID0_DEVICE_LXDAPI: // for real time stream type messages.
            switch (msgtype_subid)
            {
                case 0:
                    GetStreamData_LXDeviceAPI(uintWPARAM); // the new stream data is allocated on stStreamData
                    which is ST_STREAMDATA_LXDAPI type variable.
                    // At this point, you can use stStreamData
                    break;
            } // switch (msgtype_subid)
            break; // case MSGTYPEID0_DEVICE_LXDAPI: // for real time stream type messages.
    } // switch (msgtype_id)

    return 0;
}
```



API Functions





OpenApi

`int` OpenApi_LXDeviceAPI(`int` api_window, `int` api_selfupdate, `int` mode)

Parameter

Parameter	Available Value.	Description
<code>int</code> api_window	0 : "api window" view off. 1 : "api window" view on.(default)	api window view enable/disable.
<code>int</code> api_selfupdate	0 : "api check self update" execution. 1 ; "api check self update" no exec.(default)	api self-upadte check execution or not.
<code>int</code> mode	0 ⊕default)	Reserved

Return Values


Return	Meaning	Description
1	Success	
-1	Fail	Already succesfully called OpenApi.
-2	Fail	Windows 8/10. LXDeviceAPI fixing the os registry for USB normal communication.
-5	Fail. Duplicated call	When OpenApi called, LXDeviceAPI loading the internal used DLL from the folder LXDeviceAPI_Utility.

Code Example

```
void CLXDeviceAPI_Sample1View::OnMenuOpenapi()
{
    OpenApi_LXDeviceAPI(1,0,0);
}
```

Result

If OpenApi is called successfully, LXDeviceAPI window & tray icon appears at the bottom right.



LXDeviceAPI Window Enabled.



CloseApi

`int CloseApi_LXDeviceAPI()`

Return Values

Return Values	Meaning	Description
1	Success	
-10	Fail. Nothing to close	CloseApi should be called after OpenApi.

Code Example

```
void CLXDeviceAPI_SampleView::OnMenuCloseapi()
{
    CloseApi_LXDeviceAPI();
}
```




OpenDevice

`int` OpenDevice_LXDeviceAPI(`int` LXDeviceID, `ST_STREAMDATA_LXDAPI*` p_streamdata, `int` numsample_return, `int` mode)

Parameter

Parameter	Available Value	Description
<code>int</code> LXDeviceID	1~65535	Unique ID of the device to communicate.
<code>ST_STREAMDATA_LXDAPI*</code> p_stream_data	Address struct type of <code>ST_STREAMDATA_LXDAPI</code>	example. <code>ST_STREAMDATA_LXDAPI</code> stStreamData; <code>OpenDevice_LXDeviceAPI</code> (, &stStreamData, ,);
<code>int</code> numsample_return	1~128	The number of sampling point per one stream message.
<code>int</code> mode	0: (default)	reserved

Return Values

Return Values	Meaning	Description
>0	Success. Device Handling ID.	Application program must save the return value to call the other functions which has a parameter <code>int</code> device_handling_id
-1	Fail. Not supporting device.	
-2	Fail. Duplicated Call	Already successfully called OpenApi.
-3	Fail. No device found matching LXDeviceID.	
-5	Fail. Need to update the device firmware.	Try again after completing the device firmware update.
-10	Fail. Wrong Calling order	OpenDevice should be called after OpenApi.
-20	Fail. Wrong paramter	<code>int</code> numsample_return should be range of 1~128.




Code Example1

```
void CLXDeviceAPI_Sample1View::OnMenuOpendevice()
{
    int retv = OpenDevice_LXDeviceAPI(300, &stStreamData, NumSampleReturn, 0); // QEEG-32FX(LXDeviceID=300) Open

    if(retv > 0) // if success
    {
        m_iDeviceHandlingID = retv;
    }
    else
    {
        AfxMessageBox(_T("Fail to OpenDevice"));
    }
}
```

Result .



Device Control Panel Enabled.

Code Example2

```
/*
OpenAPI-> OpenDevice -> StartStream .
*/
void CLXDeviceAPI_Sample1View::OnMenuOneclickstart()
{
    int retv_openapi, retv_opendevide;

    retv_openapi = OpenApi_LXDeviceAPI(1, 0, 0); // API Open.

    if (retv_openapi > 0)
    {
        retv_opendevide = OpenDevice_LXDeviceAPI(300, &stStreamData, NumSampleReturn, 0); // QEEG-32FX Open
        if (retv_opendevide > 0)
        {
            m_iDeviceHandlingID = retv_opendevide;
            SetMessageDevice_LXDeviceAPI(m_iDeviceHandlingID, 0, this->m_hWnd,
WM_LXDAPI_MSGTYPEID_0, 1);

            StartStream_LXDeviceAPI(m_iDeviceHandlingID);
        }
    }
}
```



CloseDevice

```
int CloseDevice_LXDeviceAPI(int device_handling_id);
```

Parameter

Parameter	Available Value.	Description
int device_handling_id	return value of OpenDevice	

Return Values

Return Values	Meaning	Description
1	Success	
-3	Fail. No device matching device_handling_id.	
-10	Fail. Wrong calling order.	CloseDevice should be called after OpenApi.

Code Example

```
void CLXDeviceAPI_Sample1View::OnMenuClosedevice()  
{  
    CloseDevice_LXDeviceAPI(m_iDeviceHandlingID);  
}
```



StartStream

`int StartStream_LXDeviceAPI(int device_handling_id, int mode)`

Parameter

Parameter	Available Value.	Description
<code>int device_handling_id</code>	return value of OpenDevice	
<code>int mode</code>	0 : (Default)	Reserved.

Return Values

Return Values	Meaning	Description
1	Success	
-3	Fail	No device(correspond to <code>int device_handling_id</code>).
-5	Fail	Impedance measuring mode.
-6	Fail	Auto calibration mode.
-10	Fail	Wrong calling order. StartStream should be called after OpenApi.

Code Example

```
void CLXDeviceAPI_SampleView::OnMenuStartstream()
{
    StartStream_LXDeviceAPI(m_iDeviceHandlingID);
}
```



StopStream

`int StopStream_ LXDeviceAPI (int device_handling_id);`

Parameter

Parameter	Available Value.	Description
<code>int device_handling_id</code>	return value of OpenDevice	

Return Values

Return Values	Meaning	Description
1	Seccess	
-3	Fail	No device(correspond to <code>int device_handling_id</code>).
-10	Fail	Wrong calling order. StopStream should be called after OpenApi.

Code Example

```
void CLXDeviceAPI_Sample1View::OnMenuStopstream()
{
    StopStream_LXDeviceAPI(m_iDeviceHandlingID);
}
```



GetStreamData

```
int GetStreamData_LXDeviceAPI(unsigned int message_wparam);
```

Parameter

Parameter	Available Value.	Description
unsigned int message_wparam	(unsigned int) wParam.	The parameter wParam from application's message handler function.

Return Values

Return Values	Meaning	Description
1	Success	
-3	Fail	No device(correspond to int device_handling_id).
-10	Fail	Wrong calling order. GetStreamData should be called after OpenApi.

Code Example

```
afx_msg LRESULT CLXDeviceAPI_Sample1View::OnStreamData(WPARAM wParam, LPARAM lParam)
{
    unsigned int uintWPARAM = (unsigned int)wParam;

    unsigned char msgtype_id = (unsigned char)(uintWPARAM >> 8); //get the lowest 2'nd byte(message
type id).

    unsigned char msgtype_subid = (unsigned char)(uintWPARAM); //get the lowest 1st byte(message type
sub id).

    switch (msgtype_id)
    {
        case MSGTYPEID0_DEVICE_LXDAPI: // for real time stream type messages.

            switch (msgtype_subid)
            {
                case 0:

                    GetStreamData_LXDeviceAPI(uintWPARAM); // the new stream data is allocated on
stStreamData which is ST_STREAMDATA_LXDAPI type variable.

                    /// Data arrange for Wave Stream Data for plotting.
                    for (int i = 0; i < NumSampleReturn * NumChannel_Wave; i++)
                    {
                        testfloat_Wave[i] = (float)stStreamData.Wave_StreamData_CS[i];
                    }

                    /// Data arrange for Event Marking for plotting.
                    for (int idx_event = 0; idx_event < NumSampleReturn; idx_event++)
                        testfloat_Wave[NumSampleReturn * NumChannel_Wave + idx_event] =
(float)stStreamData.Event_StreamData_CS[idx_event];

                    // Test plotting the wave and event marking.
                    ACQPLOT_DLL_Array_Datain_Strip(testfloat_Wave, NumChannel_Wave + 1,
NumSampleReturn); // parameter 1 : wave data float type array. parameter 2 : total number of wave channel +
1 for Event Marking. parameter 3 : Number of Samples per one channel.

                    break;
            }
    }
}
```



```
        } // switch (msgtype_subid)

        break; // case MSGTYPEID0_DEVICE_LXDAPI: // for real time stream type messages.
    } // switch (msgtype_id)

    return 0;
}
```



EventMarkingOnStream

`int` EventMarkingOnStream_LXDeviceAPI(`int` device_handling_id, `unsigned int` event_id);

Parameter

Parameter	Available Value.	Description
<code>int</code> device_handling_id	return value of OpenDevice	
<code>unsigned int</code> event_id	1~65535.	The event_id will be marked on signal stream. When you retrieve the marked value from stream, you can recognize the what type of the event(press the key , image shown, etc.) has been marked on stream.

Return Values

Return Values	Meaning	Description
1	Success	
-3	Fail	No device(correspond to <code>int</code> device_handling_id).
-10	Fail	Wrong calling order. Should be called after OpenApi.

Code Example 1

```
// when menu clicked, let's set the event_id as 20000
void CLXDeviceAPI_Sample1View::OnMenuEventmarkingonstream()
{
    EventMarkingOnStream_LXDeviceAPI(m_iDeviceHandlingID, 20000);
}
```




Code Example 2

```
// when keyboard pressed. let's set the event_id as arrow ky down. up 30000, down 40000, left 50000, right 60000

void CLXDeviceAPI_Sample1View::OnKeyDown(UINT nChar, UINT nRepCnt, UINT nFlags)
{
    switch (nChar)
    {
        case VK_UP:
            EventMarkingOnStream_LXDeviceAPI(m_iDeviceHandlingID, 30000);
            break;

        case VK_DOWN:
            EventMarkingOnStream_LXDeviceAPI(m_iDeviceHandlingID, 40000);
            break;

        case VK_LEFT:
            EventMarkingOnStream_LXDeviceAPI(m_iDeviceHandlingID, 50000);
            break;

        case VK_RIGHT:
            EventMarkingOnStream_LXDeviceAPI(m_iDeviceHandlingID, 60000);
            break;

    }

    // TODO: Add your message handler code here and/or call default
    CView::OnKeyDown(nChar, nRepCnt, nFlags);
}
```



SendMessageDevice

```
int SendMessageDevice_LXDeviceAPI(int device_handling_id,int msgtype_id,HWND hwnd_msgrcv,int msg_id,int onoff);
```

Parameter

Parameter	Available Value.	Description
int device_handling_id	return value of OpenDevice	
int msgtype_id	MSGTYPEID0_DEVICE_LXDAPI ... MSGTYPEID9_DEVICE_LXDAPI	Predefined on LXDeviceAPI.h
HWND hwnd_msgrcv	window handle	App.'s window handle to receive the message from LXDeviceAPI.
int msg_id	WM_USER ~ WM_USER+31643 WM_APP ~ WM_APP+16383	You should define the unique message id without duplication with the other messages in your project.
int onoff	0 : message off 1 : message on	message(correspond to int msgtype_id) on/off

Return Values

Return Values	Meaning	Description
1	Success	
-3	Fail	No device(correspond to int device_handling_id).
-10	Fail	Wrong calling order. SendMessageDevice should be called after OpenApi.

Code Example

```
#define WM_STREAM_DEVICE WM_USER+203 // Define Message to get message from LXDeviceAPI.  
  
void CLXDeviceAPI_Sample1View::OnMenuSetmessagedevice()  
{  
    SendMessageDevice_LXDeviceAPI(m_iDeviceHandlingID, MSGTYPEID0_DEVICE_LXDAPI ,this->m_hwnd,  
WM_STREAM_DEVICE,1); //  
}
```



SetSampleFrequency

`int` SetSampleFrequency_LXDeviceAPI(`int` device_handling_id, `unsigned int` sample_frequency)

Parameter

Parameter	Available Value.	Description
<code>int</code> device_handling_id	return value of OpenDevice	
<code>unsigned int</code> sample_frequency	250,500,1000,2000	Sampling Frequency. unit : Hz.

Return Values

Return Values	Meaning	Description
1	Success	
-3	Fail	No device(correspond to <code>int</code> device_handling_id).
-10	Fail	Wrong calling order. SetMessageDevice should be called after OpenApi.
-20	Fail	Not supporting value of <code>unsigned int</code> sample_frequency

Code Example

```
/*
기기의 샘플링 주파수 2000Hz로 설정.
*/
void CLXDeviceAPI_Sample1View::OnSetsamplefrequency2000hz()
{
    SetSampleFrequency_LXDeviceAPI(m_iDeviceHandlingID, 2000);
}
```



SetDeviceControlPanel

`int SetDeviceControlPanel_LXDeviceAPI(int device_handling_id, int para0, int para1);`

Parameter

Parameter	Available Value.	Description
<code>int device_handling_id</code>	return value of OpenDevice	
<code>int para0</code>	0 : Select all elements.	Target elements in API's Device Control Panel.
<code>int para1</code>	0 : disable 1 : enable	Enable/disable the target elements selected by <code>int para0</code>

Return Values

Return Values	Meaning	Description
1 이상	Success	
-3	Fail	No device(correspond to <code>int device_handling_id</code>).
-10	Fail	Wrong calling order. StopStream should be called after OpenApi.

Code Example

```
void CLXDeviceAPI_Sample1View::OnMenuSetdevicecontrolpanel()
{
    static int para1=0;
    SetDeviceControlPanel_LXDeviceAPI(m_iDeviceHandlingID, 0, para1);
    // toggling para1.
    if (para1) para1 = 0;
    else para1 = 1;
}
```



GetFilterFrequency

```
int GetFilterFrequency_LXDeviceAPI(int device_handling_id,int signal_source, float * freq_hpf, float * freq_lpf, float * freq_notch);
```

Parameter

Parameter	Available Value.	Description
int device_handling_id	return value of OpenDevice	
int signal_source	0 : EEG 1 : EOG 2 : ECG	Select bio signal type to retrieve the applied filter information.
float * freq_hpf	output	High Pass Filter cut frequency. unit : Hz. The value -100 means No HPF applied.
float * freq_lpf	output	Low Pass Filter cut frequency. unit : Hz. The value -100 means no LPF applied.
float * freq_notch	output	Notch Filter frequency. unit : Hz. The value -100 means no notch filter applied.

Return Values

Return Values	Meaning	Description
1	Success	
-3	Fail	No device(correspond to int device_handling_id).
-10	Fail	Wrong calling order. GetFilterFrequency should be called after OpenApi.



Code Example


```
void CLXDeviceAPI_Sample1View::OnMenuGetfilterfrequency()
{
    float f_hpf, f_lpf, f_notch;
    CString cst_hpf, cst_lpf, cst_notch;
    // Get EEG Filter info
    if (GetFilterFrequency_LXDeviceAPI(m_iDeviceHandlingID, 0, &f_hpf, &f_lpf, &f_notch) == 1)
    {
        if (f_hpf < -1.f)        cst_hpf = _T("HPF=off, ");
        else                    cst_hpf.Format(_T("HPF=%0.1fHz, "), f_hpf);

        if (f_lpf < -1.f)        cst_lpf = _T("LPF=off, ");
        else                    cst_lpf.Format(_T("LPF=%0.1fHz, "), f_lpf);

        if (f_notch < -1.f)      cst_notch = _T("Notch=off");
        else                    cst_notch.Format(_T("Notch=%0.1fHz"), f_notch);

        AfxMessageBox(_T("EEG Filter : ") + cst_hpf + cst_lpf + cst_notch);
    }
}
```

Result





GetSampleFrequency

```
int GetSampleFrequency_LXDeviceAPI(int device_handling_id, int * sample_frequency);
```

Parameter

Parameter	Available Value.	Description
int device_handling_id	return value of OpenDevice	
int * sample_frequency	output	

Return Values

Return Values	Meaning	Description
1	Success	
-3	Fail	No device(correspond to int device_handling_id).
-10	Fail	Wrong calling order. GetFilterFrequency should be called after OpenApi.


Code Example

```
void CLXDeviceAPI_Sample1View::OnMenuGetsamplefrequency()
{
    int sample_frequency;
    CString cst;

    if (GetSampleFrequency_LXDeviceAPI(m_iDeviceHandlingID, &sample_frequency) == 1)
    {
        cst.Format(_T("Sampling Frequency = %d Hz"), sample_frequency);

        AfxMessageBox(cst);
    }
}
```

Result





GetEEGRefElectrode

`int GetEEGRefElectrode_LXDeviceAPI(int device_handling_id, int * eeg_refelectrode);`

Parameter

Parameter	Available Value.	Description
<code>int device_handling_id</code>	return value of OpenDevice	
<code>int * eeg_refelectrode</code>	0 : A1 1 : A2 2 : A1,A2 3 : Cz (ch18)	Retrieving the selected EEG reference electrode. Cz is available only if EEG Cap is introduced.

Return Values


Return Values	Meaning	Description
1	Success	
-3	Fail	No device(correspond to <code>int device_handling_id</code>).
-10	Fail	Wrong calling order. GetFilterFrequency should be called after OpenApi.

Code Example

```
void CLXDeviceAPI_Sample1View::OnMenuGeteeegrefelectrode()
{
    int eeg_refelectrode;
    CString cst_eegref;

    if (GetEEGRefElectrode_LXDeviceAPI(m_iDeviceHandlingID,&eeg_refelectrode) == 1)
    {
        if (eeg_refelectrode == 0)            cst_eegref = _T("A1");
        else if (eeg_refelectrode == 1)      cst_eegref = _T("A2");
        else if (eeg_refelectrode == 2)      cst_eegref = _T("A1+A2");
        else if (eeg_refelectrode == 3)      cst_eegref = _T("Cz(ch18)");
        AfxMessageBox(_T("EEG Reference Electode = ") + cst_eegref);
    }
}
```

Result





CheckForUpdate

`int CheckForUpdate_LXDeviceAPI(int closeifnoupdate);`

Parameter

Parameter	Available Value.	Description
<code>int closeifnoupdate</code>	0, 1	1 : Auto close update window if there is no update 0 : Don't close update window although there is no update.


Return Values

Return Values	Meaning	Description
1 이상	Success	
-10	Fail	Wrong calling order. Should be called after OpenApi.

Code Example

```
void CLXDeviceAPI_Sample1View::OnMenuCheckforupdate()
{
    CheckForUpdate_LXDeviceAPI(0); // parameter 0 : No Auto Close
    //CheckForUpdate_LXDeviceAPI(1); // parameter 1 : Auto close check program if there is no update.
}
```

Result





Appendix 1. Supporting Devices .

Device Model	LXDeviceID	Sampling Frequency (Hz)	Channel index vs. signal source
QEEG-32FX	300	250Hz, 500Hz, 1000Hz, 2000Hz	0~31 : EEG channel 1~32, 32:EOG1, 33:EOG2, 34:ECG
QEEG-64FX(8ch)	16408	250Hz, 500Hz, 1000Hz, 2000Hz	0~7 : EEG channel 1~8, 8:EOG1, 9:EOG2, 10:ECG
QEEG-64FX(16ch)	16416	250Hz, 500Hz, 1000Hz, 2000Hz	0~15 : EEG channel 1~16, 16:EOG1, 17:EOG2, 18:ECG
QEEG-64FX(24ch)	16424	250Hz, 500Hz, 1000Hz, 2000Hz	0~23 : EEG channel 1~24, 24:EOG1, 25:EOG2, 26:ECG
QEEG-64FX(32ch)	16432	250Hz, 500Hz, 1000Hz, 2000Hz	0~31 : EEG channel 1~32, 32:EOG1, 33:EOG2, 34:ECG
QEEG-64FX(40ch)	16440	250Hz, 500Hz, 1000Hz,	0~39 : EEG channel 1~40, 40:EOG1, 41:EOG2, 42:ECG
QEEG-64FX(48ch)	16448	250Hz, 500Hz, 1000Hz,	0~47 : EEG channel 1~48, 48:EOG1, 49:EOG2, 50:ECG
QEEG-64FX(56ch)	16456	250Hz, 500Hz, 1000Hz,	0~55 : EEG channel 1~56, 56:EOG1, 57:EOG2, 58:ECG
QEEG-64FX(64ch)	16464	250Hz, 500Hz, 1000Hz,	0~63 : EEG channel 1~64, 64:EOG1, 65:EOG2, 66:ECG



Revision History

Release Date	Doc. ID	Description of Change
2017-04-25	LXD64 V1	First release